



2026  
Professional  
Product Catalog

**GUM® Toothbrushes**

Good Gum Health ..... 2  
 Adult Toothbrushes..... 3-5  
 Specialty Toothbrushes ..... 6

**Children’s Oral Care Products**

Toothbrushes ..... 7-10  
 Floss Picks..... 8

**Interproximal Cleaning**

Soft-Picks®..... 11  
 Proxabrush® Go-Betweens® Cleaners..... 12  
 Proxabrush® Handles..... 13  
 Proxabrush® Refills..... 13  
 Proxabrush® Go-Betweens® Angle Cleaners . 13  
 Flossers ..... 14  
 Dental Floss ..... 15  
 Floss Accessories..... 16

**Chairside Products**

Butler® Prophy Products..... 15  
 Adult and Kids Dealer Bundle Kits..... 18  
 Imprinted Toothbrush Program..... 19  
 Oral Care Kits ..... 20  
 Orthodontic Kits ..... 20

**Accessories**

Stimulator ..... 21  
 2-in-1 Tongue Cleaner ..... 21  
 Orthodontic Wax ..... 21  
 Mouth Mirror..... 21

**SUNSTAR Diabetes Campaign**

Diabetes is a growing epidemic that needs all of our attention. Being strongly focused on helping people live healthier lives by ensuring excellence in oral health, SUNSTAR developed a toolkit to raise the awareness about the bi-directional link between oral health and diabetes and to provide guidelines for its prevention and better management. The toolkit contains separate types of information for 4 different targets: patients, oral care professionals, pharmacists and physicians. Please contact us if you would like to receive the complete toolkit.



**DID YOU KNOW**

Diabetes and periodontitis impact each other’s onset, course and outcomes. Periodontitis interacts with diabetes in ways which impact individuals with diabetes at all stages of their disease, from early glycemic control, to complications and death.

**DIABETES** increases the prevalence and severity of periodontitis, particularly when a patient has poor glycemic control.

People with diabetes have a higher prevalence of periodontitis than the general population<sup>1</sup> - Type 2 diabetes is a risk factor for periodontitis<sup>2</sup>


Individuals with diabetes are up to 3 times more likely to suffer from periodontitis than people without diabetes<sup>1</sup>

**PATIENTS WITH PERIODONTITIS & DIABETES** are impacted far worse

- Periodontitis adversely effects glycemic control and increases the risk of diabetic complication<sup>1</sup>
- Makes it harder to achieve glycemic control<sup>1</sup>
- Increases the risk of cardiovascular complications<sup>1</sup>
- Increases the mortality risk up to 8 times in people with severe periodontitis<sup>3</sup>

1. Sanz M et al, (Consensus and Guidelines IDF and EFP) Diabetes Res Clin Pract. 2017 Dec 5;  
 2. Chávarry NG et al (Meta-analysis). Oral Health Prev Dent 2009  
 3. Saremi A, Nelson RG, Tulloch-Reid M, et al. Periodontal disease and mortality in Type 2 diabetes. Diabetes Care 2005; 28(1): 27-32.

## Good Gum Health Starts with a Clean Sulcus



### Extremely Tapered Bristles (ETB)

Precision tapering enables the Extremely Tapered Bristles (ETB) to penetrate more deeply to clean up to 3.1 mm into the sulcus and clean up to 1.4 cm along the gingival margin and between the teeth.<sup>1</sup> The upper 4 mm of the bristle is tapered to a conical tip end.


Available on: *Deep Clean Technique®* and *Sensitive Clean Technique®*

#### Bristle Detail



Extremely Tapered Bristles (ETB)  
30x magnification

#### Beneath the Sulcus



3.1 mm into the sulcus<sup>1</sup>

#### Between the teeth\*\*



1.4 cm interproximal<sup>1</sup>




### Quad-Grip® Thumb Pad

The Quad-Grip® handle was designed to guide the hand to hold the brush at 45° in all quadrants so the bristles are properly positioned at the sulcus for optimal subgingival cleaning.

*The Quad-Grip® handle promotes the modified Bass technique.*

Available on: *All Technique® Toothbrushes*



### Dome Trim® Bristle Design

Dome Trim® bristles are clinically designed to clean below the gumline<sup>2,3</sup> where periodontal disease starts. The dome shape bristling structure helps promote the recommended 45° angle when using the modified Bass technique.

Available on: *Complete Care Technique®, Sensitive Care Technique®, Technique® Classic, Super Tip®, Dome Trim®, G·U·M Bright Care™*



\* DOF-0043

1. DOF-0055

2. Microbiological and Clinical Alterations from using Butler GUM toothbrushes. Reynolds,H.S.,& J.J. Zambon. J Dent Res, Vol. 76, p. 233, ID 1753 (IADR Abstracts) 1997.

3. Butler "GUM" Toothbrushes -A Microbiological and Clinical Study. Reynolds,H.S.,& J.J. Zambon. J Dent Res, Vol. 77, p. 287, ID 1452 (IADR Abstracts) 1998.

\*\* Interproximal measures from distal to mesial sides of a tooth



**TECHNIQUE® DEEP CLEAN**

- Extremely Tapered Bristles clean up to 3.1 mm into the sulcus and clean up to 1.4 cm along the gingival margin and between the teeth<sup>1</sup>
- **10x better plaque removal** below the gumline for healthier gums<sup>1</sup>

**4.3**


AVAILABLE for IMPRINT



Now with  
**Antibacterial Bristles<sup>1</sup>**



Item #	Texture	Size	Tufts
524	Soft	Full	46
525	Soft	Compact	30



**TECHNIQUE® ENAMEL CLEAN**

- Micro-feathering splits each bristle into several fine tips for **70% more enamel clean** compared to an end rounded bristle<sup>2</sup>

AVAILABLE for IMPRINT



Item #	Texture	Size	Tufts
4471	Soft	Compact	31



**TECHNIQUE® SENSITIVE CLEAN**

- For **6x better clean<sup>3</sup>** with ultra-soft bristles specially designed for sensitivity
- Extremely tapered bristles are very fine in diameter to be extra gentle on sensitive teeth and gums, while effectively cleaning below the gumline

Now with  
**Antibacterial Bristles<sup>4</sup>**

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
527	Sensitive	Compact	30



**TECHNIQUE® SENSITIVE CARE**

- The unique Quad-Grip® handle and Dome Trim® bristles, combined with ultra-delicate bristles, offer a gentle yet effective solution for plaque removal and gum care

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
516	Sensitive	Full	46
517	Sensitive	Compact	30

\* DOF-0043 1. DOF-0055 2. DOF-0024 3. DOF-0057 4. DOF-0047



## TECHNIQUE® COMPLETE CARE

- Improved cleaning in hard-to-reach areas
- Raised bristles at the tip provide greater access to posterior regions and clean hard-to-reach surfaces
- Dome Trim® bristles are clinically designed to clean below the gumline<sup>1,2</sup>

AVAILABLE for IMPRINT



Item #	Texture	Size	Tufts
590	Soft	Full	46
591	Soft	Compact	30



## DOME TRIM®

- When brushing at the dentist-recommended 45° angle, the Dome Trim® bristle design helps to effectively clean under the gumline where periodontal disease starts<sup>1,2</sup>

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
456	Soft	Full	38
457	Soft	Compact	30
459	Sensitive	Compact	30



## Summit® +

- Extremely Tapered Bristles for deeper cleaning between teeth and below the gumline

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
505	Soft	Compact	31
509	Sensitive	Compact	31



with  
Tongue  
cleaner

## MULTI-CLEAN

- Multi-level bristles with interdental cut help access hard-to-reach areas
- Tongue Cleaner to help reduce bad breath<sup>3</sup>
- Ergonomic handle with a design that improves control for a good cleaning

AVAILABLE for IMPRINT



Item #	Texture	Size	Tufts
135	Soft	Medium	46

\* DOF-0043

1. Microbiological and Clinical Alterations from using Butler GUM toothbrushes. Reynolds, H.S., & J.J. Zambon. J Dent Res, Vol. 76, p. 233, ID 1753 (IADR Abstracts) 1997.

2. Butler "GUM" Toothbrushes - A Microbiological and Clinical Study. Reynolds, H.S., & J.J. Zambon. J Dent Res, Vol. 77, p. 287, ID 1452 (IADR Abstracts) 1998.

3. Pedrazzi V, Sato S, de Mattos Mda G, Lara EH, Panzeri H. Tongue-cleaning methods. J Periodontol. 2004 Jul;75(7):1009-12.



## G·U·M BRIGHT CARE™ Toothbrush

- Special BrightLine™ Grip bristles help remove stains for **28% whiter teeth**<sup>1</sup>
- Raised Super Tip® bristles at the tip provide access to posterior regions and clean hard-to-reach surfaces
- Dome Trim® bristles are clinically designed to clean below the gumline<sup>2,3</sup>

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
4591	Soft	Compact	30

The trademark BrightLine is owned by PMM



## G·U·M BRIGHT CARE™ SENSITIVE Toothbrush

- Special BrightLine™ Grip bristles help remove stains for **28% whiter teeth**<sup>1</sup>
- Ultra Soft bristles for a gentle touch on teeth and gums
- Dome Trim® bristles are clinically designed to clean below the gumline<sup>2,3</sup>

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
4517	Sensitive	Compact	30

The trademark BrightLine is owned by PMM



## SONIC POWERED TOOTHBRUSH

- 50% more effective plaque removal and reach between teeth<sup>4</sup>
- Tested and proven to have **45x deeper access** below gum line<sup>4</sup>



Item #	Color	Texture	Size	Tufts
4100TWA	White	Soft	Compact	32
4100TBA	Black	Soft	Compact	32



## SONIC POWERED TOOTHBRUSH REFILLS

- 50% more effective plaque removal and reach between teeth<sup>4</sup>
- Tested and proven to have **45x deeper access** below gum line<sup>4</sup>




Item #	Color	Texture	Size	Tufts
4110URWA	White	Soft	Compact	32
4110URBA	Black	Soft	Compact	32

\* DOF-0043 1. DOF-0058 with the use of toothpaste.

2. Microbiological and Clinical Alterations from using Butler GUM toothbrushes. Reynolds, H.S., & J.J. Zambon. J Dent Res, Vol. 76, p. 233, ID 1753 (IADR Abstracts) 1997.

3. Butler "GUM" Toothbrushes - A Microbiological and Clinical Study. Reynolds, H.S., & J.J. Zambon. J Dent Res, Vol. 77, p. 287, ID 1452 (IADR Abstracts) 1998.

4. DOF-0011



## ORTHODONTIC Toothbrush

- This “V” trim brush facilitates cleaning around orthodontic appliances such as brackets, arch wires, buttons and ligatures. Also available as a travel brush


Item	Texture	Trim	Tufts
124	Soft	V-Trim	30



## END-TUFT Brush

- Small brush head addresses special maintenance concerns, including orthodontic bands, furcations, implants, distal of last molar, and other hard-to-reach areas


Item	Texture	Trim	Tufts
308	Soft	Tapered	7



## POST-OPERATION Toothbrush

- Ultra-soft .004 inch bristles are ideal for post surgical cleaning, gum disease, mouth irritations, extractions, implants and grafts

Item	Texture	Trim	Tufts
317	Delicate	Flat	31




153 - Regular  
125 - Orthodontic

## TRAVEL Toothbrush

- Dual-function handle is also a ventilated snap-lock cap when not in use
- Multi-level trim designs for enhanced subgingival and interproximal accesses
- Also available in an orthodontic “V” trim

Item	Texture	Trim	Tufts
153	Soft	Multi level	37
125	Soft	V-Trim	37



## DENTURE Brush

- Flat trim, firm bristle design effectively cleans denture surfaces. Tapered brush cleans smaller hard-to-reach areas
- Lever grip reduces effort for patients with limited dexterity

Item	Texture	Trim	Tufts
201	Firm	Flat	33



**AGES**  
**2+**

### Crayola PIP-SQUEAKS™

- Soft bristles are gentle on young teeth and gums
- Narrow, tapered brush head makes maneuvering in small mouths easier
- Suction cup base holds the toothbrush upright which helps keep bristles clean and reduces clutter on the counter
- Comes in four cute Pip-Squeaks™ characters



Item	Texture	Size	Tufts
232	Soft	Compact	19



**AGES**  
**3+**

### Crayola MARKER

- Raised tip bristles clean hard-to-reach back teeth
- Slim handle design provides an ergonomic hold and easier grip for small hands
- Suction cup base holds the toothbrush upright which helps to keep bristles clean and reduces clutter on the counter



Item	Texture	Size	Tufts
227	Soft	Compact	25



**AGES**  
**5+**



### Crayola TIMER LIGHT

- Flashing, one-minute timer light helps kids know when to switch rows and encourages longer brushing time
- Suction cup base holds the toothbrush upright which helps keep bristles clean and reduces clutter on the counter
- Comes in four colors with matching light
- Non-replaceable LED battery



Item	Texture	Size	Tufts
202	Soft	Compact	29



**AGES**  
**5+**

### Crayola POWER

- Oscillating action with soft bristles gently cleans teeth and massages gums
- Comes with a cap to keep brush head clean and is great for travel too!
- Suction cup base holds the toothbrush upright which helps keep bristles clean and reduces clutter on the counter
- Comes in four bold colors
- Non-replaceable battery



Item	Texture	Size	Tufts
2272	Soft	Compact	23

1. Bristles incorporate an antibacterial agent for continuous bristle protection during the recommended life of the brush. Bacterial growth that may affect the bristles is inhibited. The agent in the bristles does not protect against disease. As always, rinse your brush after use.



**AGES  
5+**



*Ultra-soft,  
tapered  
bristles*

**Crayola DEEP CLEAN**

- Ultra-soft Deep Clean bristles allow for better reach in between teeth for better cleaning<sup>1</sup>
- Thumb grip and triangular-shaped handle offer control and encourage brushing at a 45° angle
- Suction cup bottom offers easy access and upright storage



Item	Texture	Size	Tufts
234	Ultra soft	Compact	25




**Crayola Kids' Flossers**

- Shred-resistant floss with fluoride is gentle on young gums
- Ergonomic handle design makes it easy for small hands to hold
- Each pack contains all three colors and designs
- Yummy grape flavor



895 48 packs of 3



**Crayola TWISTABLES™ Flossers**

- Extra-strong, twisted floss with fluoride
- Longer handle makes it easy for parents to help kids floss
- 3 Twisted fruit flavors to choose from, to vary the fun and keep them flossing



853 Crayola™ Twistables™ Flossers 48 packs of 3

1. DOF-0006

**NEW!**

**Squad Up your growing patients with Squishmallows!**



**AGES 3+**

**Manual Toothbrush**

- Designed to be gentle and comfortable on teeth and gums, ultra-soft bristles clean thoroughly without irritation
- Tongue cleaner helps to remove bacteria from tongue gently
- Suction cup base is a fun and easy way to store your Squishmallows toothbrush



Item	Texture	Size	Tufts
4060SM1	Ultra soft	Compact	29



**AGES 5+**

**Power Toothbrush**

- Soft bristles on a small oscillating head provide a gentle, effective clean; comfortable without causing irritation
- Suction cup base keeps brush upright, clearing counters of clutter
- 2-AAA batteries included



Item	Texture	Size	Tufts
2272SM	Soft	Compact	23



**Fun brightly-colored handles and designs with popular Squishmallows friends!**



**Jazwares**

©2026 Jazwares, LLC. Original Squishmallows, Squishmallows, and all related trademarks, characters, names, trade dress, and indicia are the intellectual property of Kelly Toys Holdings, LLC.



AGES  
2+



**MONSTERZ Kids**

- Small and cushioned head adapted to the child's age protects the teeth and gums
- Soft bristles gently and effectively remove plaque
- Center bristles help emphasize a smaller amount of toothpaste for kids
- Ergonomic handle with thumb pad for stable, comfortable grip
- Suction cup base holds the toothbrush upright which helps to keep bristles clean and reduces clutter on the counter



Item	Texture	Size	Tufts
901	Ultra soft	Compact	19



AGES  
5+



**MONSTERZ Junior**


- Cushioned head adapted to the child's age protects the teeth and gums
- Soft bristles gently and effectively remove plaque
- Center bristles help emphasize a smaller amount of toothpaste for kids.
- Ergonomic handle with thumb pad for stable, comfortable grip
- Suction cup base holds the toothbrush upright which helps to keep bristles clean and reduces clutter on the counter



Item	Texture	Size	Tufts
902	Soft	Compact	29



AGES  
5+



**TECHNIQUE Kids'**

**Quad-Grip® handle teaches children the proper brushing method!**

- Cushioned Quad-Grip® handle with the Dome Trim® design helps kids place bristles at 45° so that optimal cleaning is achieved
- Comes in four different colors

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
221	Soft	Compact	24



TEEN



**G·U·M BRIGHT CARE™ Toothbrush**

**Perfect for Teens!**

- Special BrightLine™ Grip bristles help remove stains for **28% whiter teeth**<sup>1</sup>
- Raised Super Tip® bristles at the tip provide access to posterior regions and clean hard-to-reach surfaces
- Dome Trim® bristles are clinically designed to clean below the gumline<sup>2,3</sup>

AVAILABLE for IMPRINT



Item	Texture	Size	Tufts
4594	Soft	Subcompact	23

The trademark BrightLine is owned by PMM

1. DOF-0058 with the use of toothpaste  
 2. Microbiological and Clinical Alterations from using Butler GUM toothbrushes. Reynolds, H.S., & J.J. Zambon. J Dent Res, Vol. 76, p. 233, ID 1753 (IADR Abstracts) 1997.  
 3. Butler "GUM" Toothbrushes - A Microbiological and Clinical Study. Reynolds, H.S., & J.J. Zambon. J Dent Res, Vol. 77, p. 287, ID 1452 (IADR Abstracts) 1998.



*Healthy Gums  
(Non-compliant flossers)*



*Implants, Bridges  
or Crowns*



*Orthodontics*



670  
Comfort Flex  
MINT

632  
Original

634  
Wider Spaces

650  
Advanced

649  
Advanced  
TIGHT

### SOFT-PICKS® COMFORT FLEX MINT

Refreshing mint, flexible rubber bristles slide between teeth for gentle, safe and effective food & plaque removal.

### SOFT-PICKS® ORIGINAL

Convenient, comfortable, interdental cleaning. Soft flexible bristles have a tapered design to fit better in interproximal spaces to remove plaque.

### SOFT-PICKS® WIDER SPACES

25% wider than the original

### SOFT-PICKS® ADVANCED

With a more comfortable rubber grip and longer curved handle for improved maneuverability.

### SOFT-PICKS® ADVANCED TIGHT

With a 34% smaller tip diameter and 23% smaller base diameter<sup>2</sup>, the smaller size is for tighter spaces.

Clinically proven to remove plaque as effectively as string floss!

670 Soft-Picks® Comfort Flex Mint	72 packs of 4
632 Soft-Picks® Original	72 packs of 5
634 Soft-Picks® Wider Spaces	72 packs of 5
650 Soft-Picks® Advanced	110 packs of 2
649 Soft-Picks® Advanced Tight	110 packs of 2

\* DOF-0043 1. DOF-0016 2. DOF-0061. vs. Soft-Picks Advanced Regular




# INTERPROXIMAL CLEANING


# #1

Interdental Brush  
Recommended by Dental Hygienists\*


Triangular-shaped bristles remove up to 25% more plaque than conventional round bristles<sup>1</sup>



Triangular Bristles



Conventional Round Bristles



## PROXABRUSH® GO-BETWEENS® CLEANERS

- Triangular-shaped bristles remove 25% more plaque than conventional round bristles<sup>1</sup>
- Cleans posterior teeth easily with bendable neck
- Better comfort, control and grip with flexible rubberized handle
- Stays clean between uses with antibacterial bristles<sup>2</sup>
- Comes with protective cap



4.7

	◀ 2.08mm ▶ 2.08mm ▲ 0.36mm	◀ 2.30mm ▶ 2.30mm ▲ 0.48mm	◀ 2.31mm ▶ 2.79mm ▲ 0.48mm	◀ 3.70mm ▶ 4.30mm ▲ 0.67mm	◀ 3.79mm ▶ 6.86mm ▲ 0.67mm	◀ 6.10mm ▶ 6.10mm ▲ 0.84mm
Box (cello wrapped), 36/box	● Micro Tight 1312	● Ultra Tight 871	● Tight 872	● Moderate 3612	● Wide 3614	● Extra Wide 1618
Variety Pack, 36/box (6 of each individual)	●●●●●● 3615					

\* DOF-0043 1. DOF-0001

ORDER NOW


1-800-528-8537 | <http://Pro.GUMbrand.com> | [us.pcr@sunstar.com](mailto:us.pcr@sunstar.com)



# INTERPROXIMAL CLEANING

## PROXABRUSH® CLASSIC HANDLE & REFILLS

- Soft rubber grips for added comfort and control
- Triangular-shaped bristles remove 25% more plaque than conventional round bristles<sup>1</sup>
- Nylon-coated wire core provides comfort and will not scratch implants
- Wide variety of brush sizes and styles



611 Single Ended Handle (brushes not included)	12/box
605 Dual Ended Handle (brushes not included)	12/box
● 414PQF Tight Refills	18 packs of 2
● 612PQA Moderate Refills	18 packs of 2
● 614PQA Wide Refills	18 packs of 2
● 618PQA Ultra Wide Refills	18 packs of 2
● 612PQB Moderate Bulk Refills	200/box
● 614PQB Wide Bulk Refills	200/box
● 618PQB Ultra Wide Bulk Refills	200/box


## PROXABRUSH® GO-BETWEENS® ANGLE CLEANERS

- Triangular-shaped bristles remove 25% more plaque than conventional round bristles<sup>1</sup>
- Angled neck helps patients reach hard-to-clean interdental surfaces
- Tight tapered shape fits in spaces as small as 1.2 mm
- Comes with carrying case



● 3014 PROXABRUSH® GO-BETWEENS® ANGLE CLEANERS	Tight Tapered	36/box
--	---------------	--------

1. DOF-0001



*Twisted Mint®*

*Advanced Care*

*Professional Clean*

*Professional Clean PLUS*

**Twisted MINT® FLOSS PICKS**

- Twisted mint floss removes 2x more plaque<sup>2</sup> and covers 30% more surface area<sup>3</sup>
- Lightweight flexible grip uses 20% less plastic<sup>4</sup>

**Advanced CARE FLOSS PICKS MINT with Vitamin E & Fluoride**

- Durable floss coated with Vitamin E & Fluoride is designed to resist shredding or breaking<sup>5</sup>
- Lightweight flexible grip uses 20% less plastic<sup>4</sup>

**Professional CLEAN FLOSS PICKS MINT**

- Textured floss effectively removes 2x more plaque<sup>2</sup> and is designed to resist shredding or breaking
- Lightweight flexible grip uses 20% less plastic<sup>4</sup>

**Professional CLEAN PLUS FLOSS PICKS MINT**

- Specially designed rubber pick massages gums while gently removing debris
- Strong floss designed to slide easily between teeth and resist shredding or breaking

828 Twisted Mint® Floss Picks	48 packs of 3
888 Advanced Care Floss Picks Mint + Vit E & Fluoride	48 packs of 3
891 Professional Clean Floss Picks Mint	48 packs of 3
885 Professional Clean PLUS Floss Picks Mint	48 packs of 2

\* DOF-0043 1. DOF-0045 2. DOF-0049 3. DOF-0039 4. DOF-0045

**DEEP CLEAN EXPANDING<sup>®</sup>  
DENTAL FLOSS**

- Thin and easy to insert, even between tight contacts
- Floss expands during use to remove **2X more plaque**<sup>1</sup>
- Multiple filaments and texturized fibers for exceptional cleaning



2030 Deep Clean Expanding Waxed      4yd      144/box

**PRO-WEAVE<sup>®</sup>  
DENTAL FLOSS**

*(formerly known as BUTLERWEAVE<sup>®</sup>)*

- Removes 2x more plaque than the leading floss<sup>2</sup>
- Unique woven design effectively removes plaque between teeth and below the gumline
- The unique interlacing process produces a strong, smooth, shred-resistant floss



1115 PRO-WEAVE<sup>®</sup> Waxed      4yd      144/box  
 1815 PRO-WEAVE<sup>®</sup> Waxed Mint      4yd      144/box

1. DOF-0030      2. DOF-0049 Cwik J et al., Am J Dent 34(3), 2021, 123-126.

**EEZ-THRU® FLOSS THREADERS**

- Versatile nylon loop helps patients thread floss through spaces where most other products cannot reach: under bridges, between connected crowns, around orthodontic appliances, lingual wires and implants

840 Eez-Thru® Floss Threaders 1 box / 100 env. of 5



**EASYTHREAD® FLOSS**

**GUM® EasyThread® Floss is designed to make flossing around braces, bridges, and implants easy.**

- Stiff threader sections provide easy access to hard-to-clean spaces without irritating the gums
- Double-ended threader sections to maximize cleaning with each strand
- Features “puffy” floss designed for enhanced plaque removal
- Innovative, hygienic packaging reduces tangles and waste
- Office dispenser is rewindable to ensure just the right amount of floss is dispensed

3200PQ EasyThread® Floss Patient Packs 1 pkg. / 100 env. of 5  
 3210PQ EasyThread® Floss Office Dispensers 1 pkg. / 6 dispensers of 50



**FLOSSMATE FLOSS HANDLE**

- Longer, curved fingers allow access to posterior teeth and make it easy to place floss between teeth without interference from the handle
- Compatible with any floss of your choice

845 Flossmate Floss Handle 1 dozen/box




## BUTLER<sup>®</sup> Prophyiciency<sup>®</sup> Prophy Angles

The **ONLY** prophy angle that cleans and polishes without paste!

- No need to refill the cup
- No paste splatter, no paste mess
- Better visibility during cleaning
- Patient does not have to rinse and spit
- Easier cleanup between patients

Featuring a contouring cup designed to:

- Clean and polish effectively interproximally, subgingivally on all tooth surfaces with no paste required
- Clean & Polish designed to effectively clean and remove stains while being less abrasive and with no splatter<sup>1</sup>
- Polish designed for low abrasive polishing of teeth with no stains



No Paste,  
No Splatter!



Clean and Polish      Polish

1211P - Clean & Polish	200/box
1215P - Clean & Polish	50/box
1217P - Polish	200/box

## BUTLER<sup>®</sup> Eez-Touch<sup>®</sup> Disposable Prophy Angles

Soft grip, ideal comfort and control

- Soft grip provides greater control
- Reduces vibration, less fatigue
- Soft Cup and Tapered Brush



4.4



1201 Soft Cup	200/box
1203 Tapered Brush	50/box

## BUTLER<sup>®</sup> DuroPro<sup>®</sup> Disposable Prophy Angles

Ideal for your ortho patients

- Low profile
- Effective cleaning and stain removal
- Durable drive mechanism provides smooth and consistent performance



1207 Soft Cup	200/box
1205 Firm Cup	200/box
1209 Tapered Brush	50/box

1. DOF-0002



# ADULT AND KIDS DEALER BUNDLE KITS

## Enhance Your Practice with Enriched Patient Care: Premium Toothbrush Kits for Lasting Smiles.

### Customized Care



**KIT525P**  
Technique®  
Deep Clean



**KIT505P**  
Summit® +



**KIT135P**  
Multi-Clean




**KIT591P**  
Technique®  
Complete Care



**KIT4591P**  
G+U\*M  
Bright Care™



**KIT33P**  
Technique®  
Complete Care



**KIT4471P**  
Technique®  
Enamel Clean

## Empower Your Little Patients: Build Healthy Habits with GUM®

### Customized Care



**KIT35P**  
Crayola  
Marker Brush



**KIT36P**  
Crayola  
Pip-Squeaks™



**Each Kit Includes:**  
144 Toothbrushes  
144 4yd Floss tubs, or Floss Picks  
144 Printed Patient Bags

The trademark BrightLine is owned by PMM

**ORDER NOW**

1-800-528-8537 | <http://Pro.GUMbrand.com> | [us.pcr@sunstar.com](mailto:us.pcr@sunstar.com)



# IMPRINTED TOOTHBRUSHES PROGRAM

Purchase 36 Dozen GUM® Imprinted Toothbrushes, GET 36 Dozen GUM® Sample Floss or Floss Picks AND a BONUS box of any Soft-Picks® or any Proxabrush® Go-Betweens®



Patient Bags always included

## Eligible Toothbrushes

**PREMIUM**

- 525CNL - **Technique®** Deep Clean - soft/compact
- 524CNL - **Technique®** Deep Clean - soft/full
- 527CNL - **Technique®** Sensitive Clean - sensitive/compact
- 517CNL - **Technique®** Sensitive Care - sensitive/compact
- 591CNL - **Technique®** Complete Care - soft/compact
- 590CNL - **Technique®** Complete Care - soft/full
- 4471CNL - **Technique®** Enamel Clean - soft/compact
- 505CNL - **Summit®+** - soft/compact
- 509CNL - **Summit®+** - sensitive/compact

**STANDARD**

- 457CNL - **Dome Trim®** - soft/compact
- 456CNL - **Dome Trim®** - soft/full
- 465CNL - **Super Tip®** - sensitive/compact (While Supplies Last)
- 4591CNL - **G•U•M BRIGHT CARE™** - soft/compact
- 4517CNL - **G•U•M BRIGHT CARE™** - sensitive/compact
- 4594CNL - **G•U•M BRIGHT CARE™** - soft/subcompact
- 135CNL - **Multi-Clean** - soft/medium
- 221CNL - **Technique®** Kids - soft/sub-compact (KIDS)

GET FREE FLOSS OR FLOSS PICKS (minimum order 24 dozen)  
(Select one for every 12 dozen toothbrushes purchased)



GUM® Deep Clean Expanding® Floss 2030 (1 box - 144ea/box)



GUM® Pro-Weave® Floss, Waxed 1115 (1 box - 144ea/box)



GUM® Pro-Weave® Floss, Mint Waxed 1815 (1 box - 144ea/box)



GUM® Professional Clean Flossers 891 (3 boxes - 48packs/3/box)



GUM® Advanced Care Flossers 888 (3 boxes - 48 packs/3/box)



GUM® Twisted Mint® Flossers 828 (3 boxes - 48 packs/2/box)

### Use Promo Code: IMP2026

Dealer-invoiced no-charge goods will be redeemed by Sunstar. To receive your no-charge goods, send a copy of your invoice with office billing address & phone number dated between 01/01/2026 and 12/31/2026 to: us.pcr@sunstar.com or Sunstar Americas, Inc., Attn: Customer Service, 301 East Central Rd., Schaumburg, IL 60195 or fax 800-553-2014. Invoices must be received by Sunstar prior to 1/15/2027.

Scan QR Code to visit PRO GUM toothbrush page



**Orthodontic Standard Kit (124PQ)**

This kit provides the basic preventive products for every orthodontic patient.

*Kit contents:*

- GUM® Orthodontic Toothbrush
- GUM® Proxabrush® Go-Betweens® Cleaners - Tight, Moderate and Wide
- GUM® EasyThread® Floss
- GUM® Orthodontic Wax - Mint flavored



124PQ Orthodontic Standard Kit

1 dozen/box

**Orthodontic Premium Kit (ORTHKITPRB)**

This kit provides the best products an orthodontic patient can use to maintain a healthy smile.

*Kit contents:*

- GUM® Orthodontic Toothbrush
- GUM® Proxabrush® Go-Betweens® Cleaners - Wide
- GUM® Soft-Picks® Original
- GUM® Pro-Weave® Mint Waxed Floss 4 yds
- GUM® Eez-Thru® Floss Threaders
- GUM® Orthodontic Wax - Unflavored
- GUM® Toothbrush Cover
- Orthodontic Instructional Brochure
- Clear, Vented Bag with Business Card Pocket



ORTHKITPRB Orthodontic Premium Kit

1 dozen/box

**STIMULATOR**  
Handle & Refills

- Rubber tip stimulator for gingival massage and interproximal cleaning
- Slender angled neck for easy maneuverability around arches



600 Handle	1 dozen
600A Refills	3 dozen/box

**2-in-1 TONGUE CLEANER**

- Two rows of bristles reach into tiny grooves on the tongue's surface to dislodge bacteria and food particles
- Two rows of serrated scrapers gently remove bacteria and food particles while cleaning the tongue's surface
- Narrow head reduces gag reflex

Up to 50% of bad breath comes from residues on the tongue.  
75% of bad breath compounds are removed by tongue cleaning.\*




\*Pedrazzi V, et al. Tongue - Cleaning methods: a comparative clinical trial employing a toothbrush and a tongue scraper. J Periodontol 2004 Jul; 75: 1009-1012.

760	6/box
-----	-------

**ORTHODONTIC Wax**

- Clear wax adheres to orthodontic appliances to help relieve irritated tissue
- Available in unflavored (723) and mint (724) in a portable container.



723 Orthodontic Wax w/Vitamin E	2 dozen/box
724 Orthodontic Wax Mint w/Vitamin E+Aloe Vera	2 dozen/box

**MOUTH MIRROR**

- Disposable mirror with plastic handle (25mm diameter)
- Glass mirror (not foil) provides detailed reflective image



706	dozen
-----	-------



# Oral health is your mission. As your ally, we appreciate your commitment to this mission.



Every day, dental hygienists make it their job to improve each patient's health and quality of life.

We are your ally and champion on that shared mission to create overall wellness that starts in the mouth.

From continuing education to proudly being the exclusive co-sponsor of the **GUM® | RDH Award of Distinction** for the past 25 years, we will continue to find meaningful ways to support the world's Registered Dental Hygienists and the vital work they do.

**Pro.GUMbrand.com**  
**#WeAreYourAlly**

**ORDER NOW**

1-800-528-8537 | <http://Pro.GUMbrand.com> | [us.pcr@sunstar.com](mailto:us.pcr@sunstar.com)

How to Order  
**PROFESSIONAL**



**Phone Sunstar GUM**  
Call: 1-800-528-8537  
or your local GUM  
representative



**Phone Dealer**  
Call your local dealer branch  
or sales representative to place  
an order for GUM products.



**Web**  
<http://Pro.GUMbrand.com>

How to Order  
**CONSUMER**



**Call Sunstar GUM**  
1-888-777-3101  
Customer Service  
will take your order



**Web**  
[sunstarGUM.com](http://sunstarGUM.com)



**Follow us  
on Facebook**  
[@SunstarUSA](#)



**Follow us  
on Threads**  
[@SunstarUSA](#)



**Follow us  
on Instagram**  
[@SunstarUSA](#)



**Follow us  
on LinkedIn**  
[@Sunstar Americas, Inc.](#)



[Pro.GUMbrand.com](http://Pro.GUMbrand.com)

# SUNSTAR

GUM® is a registered trademark of SUNSTAR Suisse S.A.

Sunstar Americas, Inc.  
301 East Central Road, Schaumburg, IL 60195